



COMPONENTS | DEVICES | POSSIBILITIES™

High Pressure Flow Control Switch



High Pressure Flow Control Switch



SPECIFICATIONS

Part #	Description
702211001	Inline Flow Control Switch with ON/OFF Slider



WHAT IT DOES


Merit's High Pressure Flow Control Switch can be used with any application that involves liquid flow and requires inline shut off capability.

FEATURES & BENEFITS

Feature	Benefit
Audible Switch	Ensures ease of use and reliability
High Pressure	Capable of pressures up to 1200 psi in dynamic (ON) applications
On/Off Slider Switch	One-handed actuation

Sold and packaged bulk non-sterile. For sterile and other packaging options, please contact a Merit Medical OEM sales representative.

RELATED PRODUCTS



HOW TO ORDER



Inside USA +1 800.637.4839
Outside USA +1 801.208.4313
Europe +31 43 358 82 24
Or contact us at www.MeritOEM.com