



COMPONENTS | DEVICES | POSSIBILITIES™

Hemostasis Valves & Y-Adapters



FLO30™ Tuohy Borst



French Size up to 9F
Inner Lumen: 9F (.118")
Luer: Fixed

FLO30P

Access-9™



Large Bore
Hemostasis Valve
Inner Lumen: 9F (.118")
Luer: Rotating

500366002

Access-9™



Large Bore Hemostasis Valve
with 6" PVC Tubing Extension
Inner Lumen: 9F (.118")
Luer: Rotating

535001001

Access-9™



Large Bore Hemostasis Valve with
12" PVC Tubing Extension
Inner Lumen: 9F (.118")
Luer: Rotating

535003012

Access-9™



Large Bore Hemostasis Valve with
24" Polyurethane Tubing Extension
Inner Lumen: 9F (.118")
Luer: Rotating

535012024

Access-9™



Large Bore Hemostasis Valve with
12" Polyurethane Tubing Extension
Inner Lumen: 9F (.118")
Luer: Rotating

535012012

AccessPLUS™



Large Bore Hemostasis Valve
Inner Lumen: 9F (.118")
Luer: Rotating

500628002

AccessPLUS™



Large Bore Hemostasis Valve
with 12" PVC Tubing Extension
Inner Lumen: 9F (.118")
Luer: Rotating

563010001

DoublePlay™



Large Bore Double-Y Hemostasis Valve
Inner Lumen: 9F (.118")
Luer: Rotating

501075001

Passage®



Hemostasis Valve
Inner Lumen: 7.5F (.097")
Luer: Rotating

500066

Passage®



Hemostasis Valve with Slip Fit
(.179" to .209" taper)
Inner Lumen: 7.5F (.097")

K12-00457P

Honor®



Hemostasis Valve
Inner Lumen: 7F (.092")
Luer: Rotating

580195001

Honor®



Hemostasis Valve with 8" Tubing Extension
Inner Lumen: 7F (.092")
Luer: Rotating

MAP350

MBA™



Large Bore Dual Seal Hemostasis Valve
Inner Lumen: 9F (.118")

500981002

MBA™



Large Bore Hemostasis
Valve with 10" PVC
Tubing Extension

501086001


Valve Adapter



Adapter with Blue Stopcock Handle, Male with Swivel Luer
French Size: 8F or Smaller, Inner Lumen: .114"
Also available with green or red handle and cap.

501744001

SPECIFICATIONS



WHAT IT DOES

Merit Medical's hemostasis valves and Y-adapters are used to seal off and minimize fluid loss during interventional and diagnostic procedures.

FEATURES & BENEFITS

Feature	Benefit
Cap stop mechanism	Prevents cap separation from body when fully opened
From fully closed to 9 French	Accepts wide range of devices
Cap and seal engagement options	Choice of push-pull or twist
Proprietary valve technology	Leak-proof seal

CUSTOMIZATION OPTIONS

Flush port tubing lengths, connectors, bulk non-sterile, or sterile/private labeled.

RELATED PRODUCTS

Guide Wires



Catheters



HOW TO ORDER

Inside USA +1 800.637.4839

Outside USA +1 801.208.4313

Europe +31 43 358 82 24

Or contact us at www.MeritOEM.com