



COMPONENTS | DEVICES | POSSIBILITIES™

Luer Activated Valves



Part #	Description
 701000001	Male Luer Lock, PC
 702242001	Male Luer Lock, Red, PC
 702237001	Male Luer Lock, Copolyester
 702238001	Y-Port 4.0mm, PC
 702239001	Y-Port 4.0mm, Copolyester
 702236001	T-Port 6.8mm ID, PC
 702243001	Tube End 2.0mm, PC
 702243002	Tube End 2.5mm, PC
 702243003	Tube End 2.8mm, PC
 702243004	Tube End 3.7mm, PC
 701814001	Tube End 4.0mm, PC
 702240001	Bag Valve 6.8mm OD, PC
 702244001	Vial Adapter 13mm, PC
 702244002	Vial Adapter 20mm, PC
 702245001	Breakoff Valve 6.0mm, PC
 702245002	Breakoff Valve 6.6mm, PC
 702241001	Male Slip Luer, PC


WHAT IT DOES

Merit's Luer Activated Valves are designed to aspirate or inject fluids on demand. The proprietary valve facilitates multiple uses and, upon disconnection, does not require a dead-end cap to prevent leaking.

FEATURES & BENEFITS

Feature	Benefit
Proprietary Silicone Valve	Provides leak proof, swabable seal
ISO 10993 Compliance	All materials are DEHP-free and latex free
Material Bonding	Polycarbonate and copolyester materials for easy bonding

SPECIFICATIONS



RELATED PRODUCTS



HOW TO ORDER

Inside USA +1 800.637.4839
 Outside USA +1 801.208.4313
 Europe +31 43 358 82 24
 Or contact us at www.MeritOEM.com