



COMPONENTS | DEVICES | POSSIBILITIES™

# Marquis® Series Stopcocks



## 1-WAY STOPCOCKS

Catalog Number	PSI	Handle Color	Connector Type	Pictured
U1SLP	1200	Red	Slip Luer	
U1SNP	1200	Red	Swivel Nut	
U1FLP	1200	Red	Fixed Male Luer	Yes
U1RP	1200	Red	Right Rotating Male Luer	
H1FLP	1050	Blue	Fixed Male Luer	
H1RP	1050	Blue	Right Rotating Male Luer	Yes
H1SNP	1050	Blue	Swivel Nut	Yes
M1FLP	500	Light Blue	Fixed Male Luer	
M1SLP	500	Light Blue	Slip Luer	Yes
M1SNP	500	Light Blue	Swivel Nut	
S1FLP	200	White	Fixed Male Luer	Yes
S1SLP	200	White	Slip Luer	
S1SNP	200	White	Swivel Nut	

## 3-WAY STOPCOCKS

Catalog Number	PSI	Handle Color	Connector Type	Pictured
U3SLP	1200	Red	Slip Luer	
U3SNP	1200	Red	Swivel Nut	
U3FLP	1200	Red	Fixed Male Luer	
U3RRP	1200	Red	Right Rotating Male Luer	Yes
H3FLP	1050	Blue	Fixed Male Luer	
H3RLP	1050	Blue	Left Rotating Male Luer	
H3RRP	1050	Blue	Right Rotating Male Luer	
H3SLP	1050	Blue	Slip Luer	Yes
H3SNP	1050	Blue	Swivel Nut	
M3FLP	500	Light Blue	Fixed Male Luer	
M3RLP	500	Light Blue	Left Rotating Male Luer	Yes
M3RRP	500	Light Blue	Right Rotating Male Luer	
M3SLP	500	Light Blue	Slip Luer	
M3SNP	500	Light Blue	Swivel Nut	Yes
S3FP	200	White	All Female Luer	
S3FLP	200	White	Fixed Male Luer	
S3SLP	200	White	Slip Luer	
S3SNP	200	White	Swivel Nut	Yes
501641004	200	Red	Barbed Connector	
501641006	200	Green	Barbed Connector	
501641008	200	Blue	Barbed Connector	
L3SNP	50	White	Swivel Nut, Large Bore	
501855001	50	White	With Tubing Pocket .110	
501868001	50	White	With Tubing Pocket .158	Yes

## 4-WAY STOPCOCKS

Catalog Number	PSI	Handle Color	Connector Type	Pictured
U4SLP	1200	Red	Slip Luer	
U4SNP	1200	Red	Swivel Nut	
U4FLP	1200	Red	Fixed Male Luer	
U4RRP	1200	Red	Right Rotating Male Luer	
H4FLP	1050	Blue	Fixed Male Luer	
H4RRP	1050	Blue	Right Rotating Male Luer	Yes
H4SNP	1050	Blue	Swivel Nut	
M4FLP	500	Light Blue	Fixed Male Luer	
M4RRP	500	Light Blue	Right Rotating Male Luer	
M4SNP	500	Light Blue	Swivel Nut	Yes
MN4RRP	500	White	On Handle and Right Rotating	
S4FLP	200	White	Fixed Male Luer	
S4SLP	200	White	Slip Luer	
S4SNP	200	White	Swivel Nut	
L4SNP	50	White	Swivel Nut, Large Bore	

## VARIETY OF CONFIGURATIONS



Fixed Male Luer,  
1200 psi

**U1FLP**



Right Rotating Male Luer,  
1050 psi

**H1RP**



Slip Luer,  
500 psi

**M1SLP**



Fixed Male Luer,  
200 psi

**S1FLP**



Swivel Nut,  
1050 psi

**H1SNP**



Right Rotating Male Luer,  
1200 psi

**U3RRP**



Slip Luer,  
1050 psi

**H3SLP**



Left Rotating Male Luer,  
500 psi

**M3RLP**



Swivel Nut,  
200 psi

**S3SNP**



Swivel Nut,  
500 psi

**M3SNP**



With Tubing Pocket .158,  
50 psi

**501868001**



Right Rotating Male Luer,  
1050 psi

**H4RRP**



Swivel Nut,  
500 psi


**M4SNP**

## CUSTOMIZATION EXAMPLES



Contact a Merit Representative about customization options including handle color, port configuration, tubing pockets and assembly.

## SPECIFICATIONS



## WHAT IT DOES

Merit's Marquis® Series Stopcocks provide leak proof assurance during dynamic applications -- ideal for high pressure situations whether using saline or contrast.


## FEATURES & BENEFITS

Feature	Benefit
Shark Fin Handle	Unique, easy-to-grip and turn
Extended Port Length	Easy visualization of fluid & air — no interference with handles
Lumen Clarity	Easy identification of bubbles Easy to de-bubble
Variety of Configurations	Choice of directions, pressures, slip, fixed, swivel, or rotator luers
Low Torque	Easy to turn handles
Ultra High Pressure	1200 psi configurations available

## CUSTOMIZATION OPTIONS

Handle colors, custom printing, bulk non-sterile, or sterile/private labeled.

## RELATED PRODUCTS



## HOW TO ORDER



Inside USA +1 800.637.4839  
Outside USA +1 801.208.4313  
Europe +31 43 358 82 24  
Or contact us at [www.MeritOEM.com](http://www.MeritOEM.com)