




COMPONENTS | DEVICES | POSSIBILITIES™

FLO 30™ Tuohy Borst



SPECIFICATIONS

Part #	Description
FLO30P	Inner Lumen: 9F (.118") Luer: Fixed



WHAT IT DOES

Merit Medical's FLO 30™ Tuohy Borst is a single-use, proprietary valved adapter which provides a leak-proof seal during interventional and diagnostic procedures.

FEATURES & BENEFITS

Feature	Benefit
Less than 3 rotations to seal valve	Rapid engagement of seal
Cap stop mechanism	Prevents cap separation from body when fully opened
From fully closed to 9 French	Accepts wide range of devices
Proprietary valve technology	Leak-proof seal

CUSTOMIZATION OPTIONS

Bulk non-sterile or sterile/private labeled.

RELATED PRODUCTS

Guide Wires



Catheters



HOW TO ORDER

Inside USA +1 800.637.4839

Outside USA +1 801.208.4313

Europe +31 43 358 82 24

Or contact us at www.MeritOEM.com