



COMPONENTS | DEVICES | POSSIBILITIES™

VacLok®

Locking Vacuum Syringe



VacLok®

Locking Vacuum Syringe



WHAT IT DOES

Merit Medical's VacLok® is a single-use syringe designed to create a vacuum when withdrawal of a viscous material or evacuation of a device such as a balloon during a procedure is required.

FEATURES & BENEFITS

Feature	Benefit
Patented Design	Locking plunger design eliminates need to manually maintain a vacuum
Sturdy Design	Polycarbonate barrel and ABS plunger provides a stable platform
Large Wings	Easy to grip
Multiple Locking Positions	Allows for multiple vacuum and volume settings

SPECIFICATIONS

Part #	Capacity	Locking Positions	Quantity Per Case*
VAC110P	10ml	4	100
VAC120P	20ml	4	100
VAC130P	30ml	6	100
VAC160P	60ml	6	100



Plunger Material:
ABS

Barrel Material:
Polycarbonate

Plunger Tip
Material: Silicone
Luer: Fixed Male

*Non-Sterile

RELATED PRODUCTS

Inflation Devices



Stopcocks



HOW TO ORDER



Inside USA +1 800.637.4839
 Outside USA +1 801.208.4313
 Europe +31 43 358 82 24
 Or contact us at www.MeritOEM.com

CUSTOMIZING

From custom pad printing on the barrel to custom colored plungers, Merit Medical can work with your product design team to customize the VacLok to fit your needs. For sterile and other packaging options, please contact a Merit Medical OEM representative.

