



COMPONENTS | DEVICES | POSSIBILITIES™

MERIT X-CHANGE™

Syringe





WHAT IT DOES

The Merit X-Change is a single-use syringe designed to fill the open-ended syringe barrel with viscous materials or bone fillers from distal end. Alternatively, the plunger can also be removed to fill from proximal end. The open-ended barrel is threaded to accept a variety of tip configurations.


FEATURES & BENEFITS

Feature	Benefit
Sturdy Design	Reinforced wings and thumb ring provide stability when dispensing viscous material
Open-ended Barrel	Easily fill syringe
Threaded Barrel	Accepts a variety of tip configurations

Product ships bulk, non-sterile. For sterile and other packaging options, please contact a Merit Medical OEM representative.

SPECIFICATIONS

Part #	Description
EX20-MLP	20ml open-ended barrel syringe with removable fixed male luer tip
102217001	Fixed male luer tip
102218001	Dead-end cap



RELATED PRODUCTS

Inflation Devices



Stopcocks



HOW TO ORDER



Inside USA +1 800.637.4839
 Outside USA +1 801.208.4313
 Europe +31 43 358 82 24
 Or contact us at www.MeritOEM.com

CUSTOMIZING

From custom pad printing on the barrel to custom colored plungers or custom barrel tip configurations, Merit Medical can work with your product design team to customize the X-Change to fit your needs.

