

Adapters & Rotators

Custom tubing pockets and gamma stable materials are also available.

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS



Caps & Covers
Page 201



Tubing
Page 191



Luers & Connectors
Page 211

90 Degree Right Angle

Material: Polycarbonate



Clear. Pocket ID .231"-.202"
Taper

100723001



Clear. Pocket ID .172"-.152"
Taper

100723002

Male Rotating

Body Material: Polycarbonate
O-Ring Material: EPDM



Clear. Tubing Pocket ID
.180"-.205"

510104002



Gamma. Tubing Pocket ID
.136"-.139"

501178001

Male/Female and Male/Male Rotating

Body Material: Polycarbonate
O-Ring Material: EPDM



Clear. Male/Female Rotating

562013001



Clear. Male/Male Rotating

MMA100P