

Catheter Accessories

- 40 Hubs
- 41 Fixation Device
- 42 Marker Bands
- 43 Catheter Shafts

May not be available for all OEM applications.

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.


SEE RELATED PRODUCTS



Hemostasis Valves
Page 119



Sheath Introducers
Page 19



Needles
Page 11

Catheter Hubs

Material: Polycarbonate



Custom colors and printing options available.



4F Hub, White
100941200



4F Hub, Clear
100942220



5F Hub, Clear
100942221



5F Hub, White
100941205



6F Hub, White
.067" thru-hole
100941213



6F Hub, Clear
100942222



6F Hub, White
.050" thru-hole
100941211



7F Hub, White
100941215



7F Hub, Clear
100942223



8F Hub, White
100941218



8F Hub, Clear
100942224



5F Hub, White
0.045 inner lumen
Printed
590012039



6F Hub, White
0.050 inner lumen
Printed
590012057




Reverse



Reverse

Fixation Device



Revolution™ Needle-less Non-vascular
Catheter Securement Device (6F-24F)

CFG624

Marker Bands

In 2007, Merit acquired a technology from Lightek Corporation to create flexible radiopaque plastic laminate composition marker bands. This technology enables Merit to place polyurethane marker bands, impregnated with Tungsten, Platinum, or other precious metal powder, on polyurethane tubing. The marker bands are tough, flexible, and visible under fluoroscopy.

CUSTOM
BUILT TO ORDER

- Catheter shaft OD's down to 0.020"
- Wall thickness down to 0.003"
- Ideal for urethane family catheters




Platinum
Marker
Band



Used in Merit's own
Drainage Catheter Line
(see page 94)

Catheter Shafts

- Braided and non-braided
- Accepts .035 and .038 guide wires
- Standard and custom lengths
- 4, 5, 6 French sizes available



Impress®
proprietary
wire braid, also
available with
hydrophilic
coating

Performa®
stiffer shaft

Concierge™
PTFE lined,
wire braided