

Diagnostic Peripheral Catheters


May not be available for all OEM applications.

Single use. Sold packaged sterile. For non-sterile and other packaging options, please contact a Merit OEM representative.


SEE RELATED PRODUCTS



Catheter
Hubs
Page 40



Guide Wires
Page 29



Needles
Page 11

Diagnostic Peripheral Catheters

Performa®

Stiffer shaft, softer tip.

67	Amir Motarjeme®
68	Bentson
68	Berenstein
69	Cobra
70	DVS A1
70	DVS A2
71	Headhunter
72	Hockey Stick
72	Hook
72	KA
73	KA2
73	Mani
73	Mikaelsson
73	Mini Catheter
74	Modified Cerebral
74	Modified Hook Flush
74	Modified Hook Selective
75	MW2
76	Newton
77	Osborn
77	Pigtail Flush
77	Renal Double Curve
78	Reuter
78	Shepherd Flush
79	Shepherd Hook
80	Simmons
81	Straight Flush
81	Straight Selective
81	Ultra Bolus Flush
82	Ultra High Flow Pigtail
82	Vertebral
83	Vessel Sizing - Adult
84	Vessel Sizing - Pediatric


Impress®

Proprietary wire braid with integrated hub.



85	30cm Catheters
85	Amir Motarjeme
86	Bentson
86	Berenstein
87	Cobra
88	Mani
88	Headhunter
88	Hockey Stick
89	Hook
89	KA2
89	Modified Bolus Flush
90	Modified Cerebral
90	Modified Hook Flush
90	Modified Hook Selective
91	Mikaelsson
91	MultiPurpose
91	MW2
92	Newton
92	Osborn
93	Pigtail Flush
93	RBI
93	Renal Double Curve
94	Reuter
94	RIM
94	Shepherd Hook
95	Simmons
95	Straight Flush
96	Straight Selective
96	Shepherd Flush
96	Ultra Bolus Flush
97	Vertebral

Diagnostic Peripheral Catheter Shapes


- a. HEADHUNTER 1
- b. HEADHUNTER 3
- c. HEADHUNTER 5
- d. HEADHUNTER 6




- a. SIMMONS 1
- b. SIMMONS 2
- c. SIMMONS 3
- d. SIMMONS 4
- e. MODIFIED SIMMONS 1.5
- f. MODIFIED SIMMONS





- a. NEWTON 1
- b. NEWTON 2
- c. NEWTON 3
- d. NEWTON 4
- e. NEWTON 5
- f. OSBORN 1
- g. OSBORN 2



- a. BENTSON 1
- b. BENTSON 2
- c. BENTSON 3
- d. MANI





- a. BERENSTEIN
- b. VERTEBRAL
- c. MODIFIED CEREBRAL
- d. MW2
- e. MODIFIED MW2
- f. DVS A1
- g. DVS A2
- h. HOCKEY STICK





Diagnostic Peripheral Catheter Shapes



- a. KA
- b. KA2
- c. HOOK 0.8
- d. HOOK 1.0
- e. SHEPHERD HOOK 0.8
- f. SHEPHERD HOOK 1.0
- g. SHEPHERD HOOK 1.0 MOD.
- h. RENAL DOUBLE CURVE





- a. STRAIGHT SELECTIVE
- b. AMIR MOTARJEME
- c. AMIR MOTARJEME CANE
- d. AMIR MOTARJEME EXCHANGE
- e. REUTER
- f. MIKAELSSON
- g. MODIFIED HOOK FLUSH



- a. COBRA 1
- b. COBRA 2
- c. COBRA 3



- a. SHEPHERD FLUSH
- b. ULTRA BOLUS FLUSH
- c. MODIFIED BOLUS FLUSH
- d. STANDARD FLUSH
- e. STRAIGHT FLUSH
- f. MODIFIED HOOK FLUSH
- g. RBI
- h. RIM



- a. MINI CATHETER
- b. MODIFIED HOOK SELECTIVE 1
- c. MODIFIED HOOK SELECTIVE 2
- d. MODIFIED HOOK SELECTIVE 3



Performa Amir Motarjeme® Catheters

AMIR MOTARJEME

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	# of Side Ports
7747-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-
7747-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-
7747-E0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7747-F0	4	.042"/1.07mm	100cm	.035" (0.89mm)	2
5591-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-
5591-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-
5591-C3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
5591-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	2

AMIR MOTARJEME CANE

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	# of Side Ports
7747-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-
7747-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	- Bumper Tip
7747-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-
7747-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	- Bumper Tip
5591-E3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-
5591-F3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-




Catheters sold individually. Custom lengths and shapes available.

Performa Bentson Catheters


BENSTON 1

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Recommended Guide Wire
7721-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7721-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5578-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5578-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	



BENSTON 2

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Recommended Guide Wire
7721-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7721-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5578-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5578-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	



BENSTON 3

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Recommended Guide Wire
7721-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7721-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5578-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5578-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	



Performa Berenstein Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7718-C0	4	.042"/1.07mm	40cm	.035" (0.89mm)	
7718-30	4	.042"/1.07mm	40cm	.038" (0.97mm)	Bumper Tip
7718-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	Bumper Tip
7718-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	
7718-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7718-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5592-C3	5	.046"/1.17mm	40cm	.035" (0.89mm)	BERENSTEIN
5592-33	5	.046"/1.17mm	40cm	.038" (0.97mm)	
5592-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	
5592-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	
5592-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5592-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	



Catheters sold individually. Custom lengths and shapes available.

Performa Cobra Catheters

COBRA1



Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7710-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7710-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2 Bumper Tip
7710-D0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-
7710-40	4	.042"/1.07mm	65cm	.038" (0.97mm)	- Bumper Tip
5584-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
5584-D3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-
5584-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	-
5584-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	2

COBRA2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7710-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7710-E0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-
7710-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2 Bumper Tip
7710-50	4	.042"/1.07mm	65cm	.038" (0.97mm)	- Bumper Tip
5584-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
5584-E3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-
5584-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-
5584-53	5	.046"/1.17mm	65cm	.038" (0.97mm)	2

COBRA3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7710-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7710-F0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-
7710-30	4	.042"/1.07mm	65cm	.038" (0.97mm)	2 Bumper Tip
7710-60	4	.042"/1.07mm	65cm	.038" (0.97mm)	- Bumper Tip
5584-C3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
5584-F3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-
5584-33	5	.046"/1.17mm	65cm	.038" (0.97mm)	-
5584-63	5	.046"/1.17mm	65cm	.038" (0.97mm)	2



Catheters sold individually. Custom lengths and shapes available.

Performa DVS A1 Catheters (10cm Taper Tip)

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
5595-13	5	.042"/1.07mm	65cm	.035" (0.89mm)
5595-33	5	.042"/1.07mm	100cm	.035" (0.89mm)
5595-23	5	.042"/1.07mm	65cm	.038" (0.97mm)
5595-43	5	.042"/1.07mm	100cm	.038" (0.97mm)



Performa DVS A2 Catheters (15cm Taper Tip)

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
5595-A3	5	.042"/1.07mm	65cm	.035" (0.89mm)
5595-C3	5	.042"/1.07mm	100cm	.035" (0.89mm)
5595-B3	5	.042"/1.07mm	65cm	.038" (0.97mm)
5595-D3	5	.042"/1.07mm	100cm	.038" (0.97mm)



Catheters sold individually. Custom lengths and shapes available.

Performa Headhunter Catheters

HEADHUNTER 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7715-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7715-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5575-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5575-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	

HEADHUNTER 3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7715-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7715-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5575-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5575-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	

HEADHUNTER 5

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7715-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7715-50	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5575-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5575-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	

HEADHUNTER 6

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7715-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7715-60	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5575-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5575-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	



Catheters sold individually. Custom lengths and shapes available.

Performa Hockey Stick Catheters

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Recommended Guide Wire	
7752-10	4	.042"/1.07mm	100cm	4.5cm	.038" (0.97mm)	Bumper Tip
7543-13	5	.046"/1.17mm	100cm	4.5cm	.038" (0.97mm)	
7541-13	6	.054"/1.37mm	100cm	4.5cm	.038" (0.97mm)	



Performa Hook Catheters

HOOK 0.8

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7711-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	Bumper Tip
7711-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	
5585-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5585-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	

HOOK 1.0

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7711-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	Bumper Tip
7711-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	
5585-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5585-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



Performa KA Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
5582-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	High Flow Non-braided
5582-33	5	.052"/1.32mm	65cm	.038" (0.97mm)	



Catheters sold individually. Custom lengths and shapes available.

Performa KA2 Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
5582-13	5	.046"/1.17mm	40cm	.038" (0.97mm)
5582-43	5	.046"/1.17mm	65cm	.038" (0.97mm)
7514-13	6	.054"/1.37mm	40cm	.038" (0.97mm)
7514-23	6	.054"/1.37mm	65cm	.038" (0.97mm)



Performa Mani Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7716-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7716-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5579-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5579-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	



Performa Mikaelsson Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7719-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	
7719-10	4	.042"/1.07mm	80cm	.038" (0.97mm)	Bumper Tip
5594-13	5	.046"/1.17mm	80cm	.038" (0.97mm)	



Performa Mini Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
42035-MINI	4	.041"/1.04mm	20cm	.035" (0.89mm)	Non-braided



Catheters sold individually. Custom lengths and shapes available.

Performa Modified Cerebral Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
7723-A0	4	.042"/1.07mm	90cm	.035" (0.89mm)
7723-10	4	.042"/1.07mm	90cm	.038" (0.97mm)
5588-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)
5588-13	5	.046"/1.17mm	100cm	.038" (0.97mm)

Bumper Tip



Performa Modified Hook Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	Port Type
7603-A1	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	Non-braided
7603-B1	5	.052"/1.32mm	90cm	.035" (0.89mm)	6	Non-braided
7603-C1	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	Non-braided
7603-D1	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	Non-braided



Performa Modified Hook Selective Catheters

MODIFIED HOOK 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
5585-C3	5	.046"/1.17mm	80cm	.035" (0.89mm)



MODIFIED HOOK 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
5585-D3	5	.046"/1.17mm	80cm	.035" (0.89mm)



MODIFIED HOOK 3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
5585-E3	5	.046"/1.17mm	80cm	.035" (0.89mm)




Catheters sold individually. Custom lengths and shapes available.

Performa MW2 Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7745-10	4	.042"/1.07mm	40cm	.038" (0.97mm)	2	Bumper Tip
5590-13	5	.046"/1.17mm	40cm	.038" (0.97mm)		
5590-23	5	.046"/1.17mm	65cm	.038" (0.97mm)		
5590-33	5	.046"/1.17mm	100cm	.038" (0.97mm)		

MODIFIED MW2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7745-A0	4	.042"/1.07mm	40cm	.035" (0.89mm)	2
7745-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7745-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	2




Catheters sold individually. Custom lengths and shapes available.

Performa Newton Catheters

NEWTON 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7722-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7722-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5577-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	Bumper Tip
5577-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	



NEWTON 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7722-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7722-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5577-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	Bumper Tip
5577-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	




NEWTON 3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7722-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7722-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5577-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	Bumper Tip
5577-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	




NEWTON 4

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7722-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7722-40	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5577-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	Bumper Tip
5577-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	



NEWTON 5

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7722-E0	4	.042"/1.07mm	100cm	.035" (0.89mm)	Bumper Tip
7722-50	4	.042"/1.07mm	100cm	.038" (0.97mm)	
5577-E3	5	.046"/1.17mm	100cm	.035" (0.89mm)	Bumper Tip
5577-53	5	.046"/1.17mm	100cm	.038" (0.97mm)	



Catheters sold individually. Custom lengths and shapes available.

Performa Osborn Catheters

OSBORN 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7746-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7746-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip

OSBORN 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7746-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	
7746-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	Bumper Tip
7746-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7746-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip




Performa Pigtail Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7749-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	10	Braided
7749-40	4	.042"/1.07mm	80cm	.038" (0.97mm)	10	Braided



Performa Renal Double Curve Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7713-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7713-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7713-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	Bumper Tip
7713-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	Bumper Tip
5587-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5587-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5587-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



Catheters sold individually. Custom lengths and shapes available.

Performa Reuter Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7714-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	2	
7714-10	4	.042"/1.07mm	80cm	.038" (0.97mm)	2	Bumper Tip
7714-B0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7714-20	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	Bumper Tip
5589-A3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5589-B3	5	.046"/1.17mm	80cm	.035" (0.89mm)	2	
5589-13	5	.046"/1.17mm	80cm	.038" (0.97mm)	-	
5589-23	5	.046"/1.17mm	80cm	.038" (0.97mm)	2	



Performa Shepherd Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7602-A1	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	Non-braided
7602-B1	5	.052"/1.32mm	90cm	.035" (0.89mm)	6	Non-braided
7602-C1	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	Non-braided
7602-D1	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	Non-braided



Catheters sold individually. Custom lengths and shapes available.

Performa Shepherd Hook Catheters

SHEPHERD HOOK 0.8


Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7712-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7712-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2 Bumper Tip
5586-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
5586-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	2

SHEPHERD HOOK 1.0

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7712-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7712-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2 Bumper Tip
5586-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
5586-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2

SHEPHERD HOOK 1.0 MODIFIED

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
7712-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2
7712-D0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-
7712-30	4	.042"/1.07mm	65cm	.038" (0.97mm)	2 Bumper Tip
7712-40	4	.042"/1.07mm	65cm	.038" (0.97mm)	- Bumper Tip
5586-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	2



Catheters sold individually. Custom lengths and shapes available.

Performa Simmons Catheters

SIMMONS 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7720-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7720-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5576-F3	5	.046"/1.17mm	65cm	.035" (0.89mm)	
5576-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5576-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	



SIMMONS 1

SIMMONS 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7720-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7720-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5576-G3	5	.046"/1.17mm	65cm	.035" (0.89mm)	
5576-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5576-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	



SIMMONS 2

SIMMONS 3


Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7720-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7720-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5576-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5576-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	



SIMMONS 3

SIMMONS 4

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7720-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7720-40	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5576-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5576-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	



SIMMONS 4

MODIFIED SIMMONS 1.5

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
7720-F0	4	.042"/1.07mm	100cm	.035" (0.89mm)
7720-E0	4	.042"/1.07mm	80cm	.035" (0.89mm)
5576-E3	5	.046"/1.17mm	100cm	.035" (0.89mm)



MODIFIED SIMMONS 1.5

MODIFIED SIMMONS

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7720-53	4	.042"/1.07mm	80cm	.038" (0.97mm)	Bumper Tip



MODIFIED SIMMONS

Catheters sold individually. Custom lengths and shapes available.

Performa Straight Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7673-11	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	Non-braided
7673-21	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	Non-braided
7673-31	5	.052"/1.32mm	65cm	.035" (0.89mm)		Non-braided
7673-41	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	Non-braided
7673-51	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	Non-braided
7673-61	5	.052"/1.32mm	65cm	.038" (0.97mm)		Non-braided

STRAIGHT



Performa Straight Selective Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7748-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7748-13	4	.042"/1.07mm	110cm	.038" (0.97mm)	Bumper Tip

STRAIGHT



Performa Ultra Bolus Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7674-11	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	Non-braided
7674-21	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	Non-braided
7674-31	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	Non-braided
7674-41	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	Non-braided

ULTRA BOLUS FLUSH



Catheters sold individually. Custom lengths and shapes available.

Performa Ultra High Flow Pigtail Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
7602-41	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	Non-braided
7602-51	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	Non-braided
7602-61	5	.052"/1.32mm	110cm	.035" (0.89mm)	10	Non-braided
7602-71	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	Non-braided
7602-81	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	Non-braided
7602-91	5	.052"/1.32mm	110cm	.038" (0.97mm)	10	Non-braided

Performa Vertebral Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
7717-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	
7717-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	Bumper Tip
7717-C0	4	.042"/1.07mm	80cm	.035" (0.89mm)	
7717-30	4	.042"/1.07mm	80cm	.038" (0.97mm)	Bumper Tip
7717-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	
7717-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	Bumper Tip
5583-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	
5583-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	
5583-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	
5583-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	



Catheters sold individually. Custom lengths and shapes available.

Performa Vessel Sizing Catheters - Adult

All side ports between tip and first band

2 BAND



Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	Bands and Spacing	Style
7742-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Pigtail
7742-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Pigtail
7744-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Straight
7744-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Straight
7743-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	UltraBolus
7743-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Ultra Bolus
7755-2M65	4	.041"/1.04mm	65cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook
7755-2M100	4	.041"/1.04mm	100cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook
7602-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
7602-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
7674-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	UltraBolus
7674-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	UltraBolus
7673-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
7673-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
7603-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook
7603-2M100	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook

11 BAND

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	Bands and Spacing	Style
7602-11M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	11 (1cm apart)	Pigtail
7602-11M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	11 (1cm apart)	Pigtail
7673-11M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	11 (1cm apart)*	Straight
7673-11M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	11 (1cm apart)*	Straight
7603-11M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	11 (1cm apart)	Mod. Hook
7603-11M100	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	11 (1cm apart)	Mod. Hook

20 BAND

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	Bands and Spacing	Style
7742-20M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	20 (1cm apart)*	Pigtail
7742-20M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	20 (1cm apart)*	Pigtail
7602-20M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
7602-20M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
7602-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	5	20 (1cm apart)*	Pigtail
7602-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	5	20 (1cm apart)*	Pigtail
7673-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 (1cm apart)*	Straight
7673-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 (1cm apart)*	Straight
7603-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	20 (1cm apart)*	Mod. Hook
7603-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	20 (1cm apart)*	Mod. Hook



Catheters sold individually. Custom lengths and shapes available. *"SH" suffix identifies marker bands that are 16mm closer to the pigtail curve.

Performa Vessel Sizing Catheters - Pediatric

All side ports between bands

2 BAND

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	Bands and Spacing	Style
56535AVS2PIG	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
510035AVS2PIG	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
56535AVS2STR	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
510035AVS2STR	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
56535AVS2MHK	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Mod Hook
510035AVS2MHK	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Mod Hook

20 BAND


Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	Bands and Spacing	Style
46535AVS20MHK	4	041"/1.04mm	65cm	035" (0.89mm)	6	20 (1cm apart)	Mod Hook
56535AVS20PIG	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
510035AVS20PIG	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
56535AVS20MHK	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook
510035AVS20MHK	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook



Catheters sold individually. Custom lengths and shapes available.

Impress 30cm Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	of Side Ports
43035CB1	4	.040"/1.02mm	30cm	.035" (0.89mm)	
43035BER	4	.040"/1.02mm	30cm	.035" (0.89mm)	
43035KA2	4	.040"/1.02mm	30cm	.035" (0.89mm)	
43035STR	4	.040"/1.02mm	30cm	.035" (0.89mm)	
43035RIM	4	.040"/1.02mm	30cm	.035" (0.89mm)	
4303510PIG-NB	4	.040"/1.02mm	30cm	.035" (0.89mm)	10
53035CB1	5	.046"/1.17mm	30cm	.035" (0.89mm)	
53035BER	5	.046"/1.17mm	30cm	.035" (0.89mm)	
53035KA2	5	.046"/1.17mm	30cm	.035" (0.89mm)	
53035STR	5	.046"/1.17mm	30cm	.035" (0.89mm)	
5303510PIG-NB	5	.040"/1.02mm	30cm	.035" (0.89mm)	10



Impress Amir Motarjeme® Catheters

AMIR MOTARJEME

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46535MOT	4	.040"/1.02mm	65cm	.035" (0.89mm)
410035MOT	4	.040"/1.02mm	100cm	.035" (0.89mm)
56535MOT	5	.046"/1.17mm	65cm	.035" (0.89mm)
510035MOT	5	.046"/1.17mm	100cm	.035" (0.89mm)

AMIR MOTARJEME CANE

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46535MOTC	4	.040"/1.02mm	65cm	.035" (0.89mm)
410035MOTC	4	.040"/1.02mm	100cm	.035" (0.89mm)




Impress Bentson Catheters

BENTSON 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46535B1	4	.040"/1.02mm	65cm	.035" (0.89mm)
410035B1	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038B1	4	.040"/1.02mm	100cm	.038" (0.97mm)
56535B1	5	.046"/1.17mm	65cm	.035" (0.89mm)
510035B1	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038B1	5	.046"/1.17mm	100cm	.038" (0.97mm)
410038B1-H	4	-	100cm	.038" (0.97mm)
510038B1-H	5	-	100cm	.038" (0.97mm)

BENTSON 2



40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating

NEW

BENTSON 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410035B2	4	.040"/1.02mm	100cm	.035" (0.89mm)
510035B2	5	.046"/1.17mm	100cm	.035" (0.89mm)
410038B2-H	4	-	100cm	.038" (0.97mm)
510038B2-H	5	-	100cm	.038" (0.97mm)

BENTSON 1



40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating

NEW

Impress Berenstein Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46535BER	4	.040"/1.02mm	65cm	.035" (0.89mm)
46538BER	4	.040"/1.02mm	65cm	.038" (0.97mm)
410035BER	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038BER	4	.040"/1.02mm	100cm	.038" (0.97mm)
412538BER	4	.040"/1.02mm	125cm	.038" (0.97mm)
56535BER	5	.046"/1.17mm	65cm	.035" (0.89mm)
56538BER	5	.046"/1.17mm	65cm	.038" (0.97mm)
510035BER	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038BER	5	.046"/1.17mm	100cm	.038" (0.97mm)
512538BER	5	.046"/1.17mm	125cm	.038" (0.97mm)
46538BER-H	4	-	65cm	.038" (0.97mm)
410038BER-H	4	-	100cm	.038" (0.97mm)
56538BER-H	5	-	65cm	.038" (0.97mm)
510038BER-H	5	-	100cm	.038" (0.97mm)

25cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating
25cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating

NEW

BERENSTEIN



Catheters sold individually. Custom lengths and shapes available.

Impress Cobra Catheters

COBRA 1 (side holes)

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
465352CB1	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
465382CB1	4	.040"/1.02mm	65cm	.038" (0.97mm)	2
565352CB1	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
565382CB1	5	.046"/1.17mm	65cm	.038" (0.97mm)	2

COBRA 1 (without side holes)

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
46538CB1	4	.040"/1.02mm	65cm	.038" (0.97mm)	
410038CB1	4	.040"/1.02mm	100cm	.038" (0.97mm)	
56538CB1	5	.046"/1.17mm	65cm	.038" (0.97mm)	
510038CB1	5	.046"/1.17mm	100cm	.038" (0.97mm)	
46538CB1-H	4	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
56538CB1-H	5	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating




COBRA 2 (side holes)

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
465352CB2	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
465382CB2	4	.040"/1.02mm	65cm	.038" (0.97mm)	2
565352CB2	5	.046"/1.17mm	65cm	.035" (0.89mm)	2
565382CB2	5	.046"/1.17mm	65cm	.038" (0.97mm)	2

COBRA 2 (without side holes)

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
46535CB2	4	.040"/1.02mm	65cm	.035" (0.89mm)	
46538CB2	4	.040"/1.02mm	65cm	.038" (0.97mm)	
410038CB2	4	.040"/1.02mm	100cm	.038" (0.97mm)	
56538CB2	5	.046"/1.17mm	65cm	.038" (0.97mm)	
510038CB2	5	.046"/1.17mm	100cm	.038" (0.97mm)	
46538CB2-H	4	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
410038CB2-H	4	-	100cm	.038" (0.97mm)	40cm of Legato™ Hydrophilic Coating
56538CB2-H	5	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
510038CB2-H	5	-	100cm	.038" (0.97mm)	40cm of Legato™ Hydrophilic Coating



Catheters sold individually. Custom lengths and shapes available.

Impress Mani Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410035MANI	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038MANI	4	.040"/1.02mm	100cm	.038" (0.97mm)
510035MANI	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038MANI	5	.046"/1.17mm	100cm	.038" (0.97mm)
410038MANI-H	4	-	100cm	.038" (0.97mm)
510038MANI-H	5	-	100cm	.038" (0.97mm)



MANI

40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating



Impress Headhunter Catheters

HEADHUNTER 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410035HH1	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038HH1	4	.040"/1.02mm	100cm	.038" (0.97mm)
412535HH1	4	.040"/1.02mm	125cm	.035" (0.89mm)
510035HH1	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038HH1	5	.046"/1.17mm	100cm	.038" (0.97mm)
512535HH1	5	.046"/1.17mm	125cm	.035" (0.89mm)
410038HH1-H	4	-	100cm	.038" (0.97mm)
412538HH1-H	4	-	125cm	.038" (0.97mm)
510038HH1-H	5	-	100cm	.038" (0.97mm)
512538HH1-H	5	-	125cm	.038" (0.97mm)



HEADHUNTER 1

40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating



HEADHUNTER 3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410035HH3	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038HH3	4	.040"/1.02mm	100cm	.038" (0.97mm)
510035HH3	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038HH3	5	.046"/1.17mm	100cm	.038" (0.97mm)



HEADHUNTER 3

Impress Hockey Stick Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410038HS	4	.040"/1.02mm	100cm	4.5cm .038" (0.97mm)
510038HS	5	.046"/1.17mm	100cm	4.5cm .038" (0.97mm)



HOCKEY STICK

Catheters sold individually. Custom lengths and shapes available.

Impress Hook Catheters

HOOK 0.8


Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
465352HK0-8	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
565352HK0-8	5	.046"/1.17mm	65cm	.035" (0.89mm)	2

HOOK 1.0

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
465352HK1-0	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
565352HK1-0	5	.046"/1.17mm	65cm	.035" (0.89mm)	2




Impress KA2 Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
44038KA2	4	.040"/1.02mm	40cm	.038" (0.97mm)	 KA2 25cm of Legato™ Hydrophilic Coating
46538KA2	4	.040"/1.02mm	65cm	.038" (0.97mm)	
54038KA2	5	.046"/1.17mm	40cm	.038" (0.97mm)	
56538KA2	5	.046"/1.17mm	65cm	.038" (0.97mm)	
44038KA2-H	4	-	40cm	.038" (0.97mm)	
46538KA2-H	4	-	65cm	.038" (0.97mm)	
54038KA2-H	5	-	40cm	.038" (0.97mm)	
56538KA2-H	5	-	65cm	.038" (0.97mm)	

NEW

Impress Modified Bolus Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
4653810MUBF-NB	4	.040"/1.02mm	65cm	.038" (0.97mm)	10	Non-braided
5653810MUBF-NB	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	Non-braided



Catheters sold individually. Custom lengths and shapes available.

Impress Modified Cerebral Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410035MC	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038MC	4	.040"/1.02mm	100cm	.038" (0.97mm)
510035MC	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038MC	5	.046"/1.17mm	100cm	.038" (0.97mm)

MODIFIED CEREBRAL



Impress Modified Hook Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
465356MHK-NB	4	.040"/1.02mm	65cm	.035" (0.89mm)	6	Non-braided
490356MHK-NB	4	.040"/1.02mm	90cm	.035" (0.89mm)	6	Non-braided
565356MHK-NB	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	Non-braided
565386MHK-NB	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	Non-braided
590386MHK-NB	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	Non-braided



MODIFIED HOOK FLUSH

Impress Modified Hook Selective Catheters

MODIFIED HOOK 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
48035MHK1	4	.040"/1.02mm	80cm	.035" (0.89mm)
58035MHK1	5	.046"/1.17mm	80cm	.035" (0.89mm)



MODIFIED HOOK 1

MODIFIED HOOK 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
48035MHK2	4	.040"/1.02mm	80cm	.035" (0.89mm)
58035MHK2	5	.046"/1.17mm	80cm	.035" (0.89mm)



MODIFIED HOOK 2

MODIFIED HOOK 3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
48035MHK3	4	.040"/1.02mm	80cm	.035" (0.89mm)
58035MHK3	5	.046"/1.17mm	80cm	.035" (0.89mm)



MODIFIED HOOK 3

Catheters sold individually. Custom lengths and shapes available.

Impress Mikaelsson Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
48038MIK	4	.040"/1.02mm	80cm	.038" (0.97mm)
58038MIK	5	.046"/1.17mm	80cm	.038" (0.97mm)



Impress Multipurpose A1 Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410038MPA1	4	.040"/1.02mm	100cm	.038" (0.97mm)
412538MPA1	4	.040"/1.02mm	125cm	.038" (0.97mm)
510038MPA1	5	.046"/1.17mm	100cm	.038" (0.97mm)
512538MPA1	5	.046"/1.17mm	125cm	.038" (0.97mm)
410038MPA1-H	4	-	100cm	.038" (0.97mm)
510038MPA1-H	5	-	100cm	.038" (0.97mm)

40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating



Impress MW2 Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
44038MW2	4	.040"/1.02mm	40cm	.038" (0.97mm)
56538MW2	5	.046"/1.17mm	65cm	.038" (0.97mm)



Catheters sold individually. Custom lengths and shapes available.

Impress Newton Catheters

NEWTON 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410038N1	4	.040"/1.02mm	100cm	.038" (0.97mm)
510038N1	5	.046"/1.17mm	100cm	.038" (0.97mm)



NEWTON 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410038N2	4	.040"/1.02mm	100cm	.038" (0.97mm)
510038N2	5	.046"/1.17mm	100cm	.038" (0.97mm)



NEWTON 3

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410038N3	4	.040"/1.02mm	100cm	.038" (0.97mm)
510038N3	5	.046"/1.17mm	100cm	.038" (0.97mm)



NEWTON 4

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410038N4	4	.040"/1.02mm	100cm	.038" (0.97mm)
510038N4	5	.046"/1.17mm	100cm	.038" (0.97mm)



Impress Osborn Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
410035OSB1	4	.040"/1.02mm	100cm	.035" (0.89mm)
510035OSB1	5	.046"/1.17mm	100cm	.035" (0.89mm)



Catheters sold individually. Custom lengths and shapes available.

Impress Pigtail Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
4303510PIG-NB	4	.040"/1.02mm	30cm	.035" (0.89mm)	10	Non-braided
4653510PIG-NB	4	.040"/1.02mm	65cm	.035" (0.89mm)	10	Non-braided
4903510PIG-NB	4	.040"/1.02mm	90cm	.035" (0.89mm)	10	Non-braided
41103510PIG-NB	4	.040"/1.02mm	110cm	.035" (0.89mm)	10	Non-braided
5303510PIG-NB	5	.052"/1.32mm	30cm	.035" (0.89mm)	10	Non-braided
5653510PIG-NB	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	Non-braided
5903510PIG-NB	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	Non-braided
51103510PIG-NB	5	.052"/1.32mm	110cm	.035" (0.89mm)	10	Non-braided
5653810PIG-NB	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	Non-braided
5903810PIG-NB	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	Non-braided
51103810PIG-NB	5	.052"/1.32mm	110cm	.038" (0.97mm)	10	Non-braided




Impress RBI (renal, bifurcation, iliac) Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
465354RBI	4	.040"/1.02mm	65cm	.035" (0.89mm)	4	
565354RBI	5	.046"/1.17mm	65cm	.035" (0.89mm)	4	
56535RBI	5	.046"/1.17mm	65cm	.035" (0.89mm)		
510035RBI	5	.046"/1.17mm	100cm	.035" (0.89mm)		
56538RBI-H	5	-	65cm	.038" (0.97mm)		25cm of Legato™ Hydrophilic Coating



Impress Renal Double Curve Catheters


Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
46535RC	4	.040"/1.02mm	65cm	.035" (0.89mm)	
465352RC	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
565352RC	5	.046"/1.17mm	65cm	.035" (0.89mm)	2



Catheters sold individually. Custom lengths and shapes available.

Impress Reuter Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
48038RE	4	.040"/1.02mm	80cm	.038" (0.97mm)
58035RE	5	.046"/1.17mm	80cm	.035" (0.89mm)
58038RE	5	.046"/1.17mm	80cm	.038" (0.97mm)



REUTER

Impress RIM Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46535RIM	4	.040"/1.02mm	65cm	.035" (0.89mm)
56535RIM	5	.046"/1.17mm	65cm	.035" (0.89mm)



RIM

Impress Shepherd Hook Catheters

SHEPHERD HOOK 0.8

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
465352SH0-8	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
565352SH0-8	5	.046"/1.17mm	65cm	.035" (0.89mm)	2

SHEPHERD HOOK 1.0

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports
465352SH1-0	4	.040"/1.02mm	65cm	.035" (0.89mm)	2
565352SH1-0	5	.046"/1.17mm	65cm	.035" (0.89mm)	2



SHEPHERD HOOK 0.8




SHEPHERD HOOK 1.0

Catheters sold individually. Custom lengths and shapes available.

Impress Simmons Catheters

SIMMONS 1

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46535SIM1	4	.040"/1.02mm	65cm	.038" (0.97mm)
46538SIM1	4	.040"/1.02mm	65cm	.038" (0.97mm)
410035SIM1	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038SIM1	4	.040"/1.02mm	100cm	.038" (0.97mm)
56535SIM1	5	.040"/1.02mm	65cm	.038" (0.97mm)
56538SIM1	5	.046"/1.17mm	65cm	.038" (0.97mm)
510035SIM1	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038SIM1	5	.046"/1.17mm	100cm	.038" (0.97mm)
410038SIM1-H	4	-	100cm	.038" (0.97mm)
510038SIM1-H	5	-	100cm	.038" (0.97mm)



SIMMONS 1

40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating



SIMMONS 2

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
46538SIM2	4	.040"/1.02mm	65cm	.038" (0.97mm)
410035SIM2	4	.040"/1.02mm	100cm	.035" (0.89mm)
410038SIM2	4	.040"/1.02mm	100cm	.038" (0.97mm)
56538SIM2	5	.046"/1.17mm	65cm	.038" (0.97mm)
510035SIM2	5	.046"/1.17mm	100cm	.035" (0.89mm)
510038SIM2	5	.046"/1.17mm	100cm	.038" (0.97mm)
410038SIM2-H	4	-	100cm	.038" (0.97mm)
510038SIM2-H	5	-	100cm	.038" (0.97mm)




SIMMONS 2

40cm of Legato™ Hydrophilic Coating
40cm of Legato™ Hydrophilic Coating



MODIFIED SIMMONS

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire
48035MODS	4	.040"/1.02mm	80cm	.035" (0.89mm)
48038MODS	4	.040"/1.02mm	80cm	.038" (0.97mm)



MODIFIED SIMMONS

Impress Straight Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
4653510STF-NB	4	.040"/1.02mm	65cm	.035" (0.89mm)	10	Non-braided
4903510STF-NB	4	.040"/1.02mm	90cm	.035" (0.89mm)	10	Non-braided
5653510STF-NB	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	Non-braided
5903510STF-NB	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	Non-braided

STRAIGHT



Catheters sold individually. Custom lengths and shapes available.

Impress Straight Selective Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
46538STS	4	.040"/1.02mm	65cm	.038" (0.97mm)	
411035STS	4	.040"/1.02mm	110cm	.035" (0.89mm)	
411038STS	4	.040"/1.02mm	110cm	.038" (0.97mm)	
56538STS	5	.046"/1.17mm	65cm	.038" (0.97mm)	
511035STS	5	.046"/1.17mm	110cm	.035" (0.89mm)	
511038STS	5	.046"/1.17mm	110cm	.038" (0.97mm)	
46538STS-H	4	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
410038STS-H	4	-	100cm	.038" (0.97mm)	40cm of Legato™ Hydrophilic Coating
56538STS-H	5	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
510038STS-H	5	-	100cm	.038" (0.97mm)	40cm of Legato™ Hydrophilic Coating

STRAIGHT SELECTIVE



Impress Shepherd Flush Catheters


Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
465356SHP-NB	4	.040"/1.02mm	65cm	.035" (0.89mm)	6	Non-braided
490356SHP-NB	4	.040"/1.02mm	90cm	.035" (0.89mm)	6	Non-braided
565356SHP-NB	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	Non-braided



SHEPHERD FLUSH

Impress Ultra Bolus Flush Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# of Side Ports	
4653510UBF-NB	4	.040"/1.02mm	65cm	.035" (0.89mm)	10	Non-braided
4903510UBF-NB	4	.040"/1.02mm	90cm	.035" (0.89mm)	10	Non-braided
5653510UBF-NB	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	Non-braided
5903510UBF-NB	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	Non-braided



ULTRA BOLUS FLUSH

Catheters sold individually. Custom lengths and shapes available.

Impress Vertebral Catheters

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	
46538VER	4	.040"/1.02mm	65cm	.038" (0.97mm)	
410035VER	4	.040"/1.02mm	100cm	.035" (0.89mm)	
410038VER	4	.040"/1.02mm	100cm	.038" (0.97mm)	
412535VER	4	.040"/1.02mm	125cm	.035" (0.89mm)	
56538VER	5	.046"/1.17mm	65cm	.038" (0.97mm)	
510035VER	5	.046"/1.17mm	100cm	.035" (0.89mm)	
510038VER	5	.046"/1.17mm	100cm	.038" (0.97mm)	
512535VER	5	.046"/1.17mm	125cm	.035" (0.89mm)	
56538VER-H	5	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
510038VER-H	5	-	100cm	.038" (0.97mm)	40cm of Legato™ Hydrophilic Coating
46538VER-H	4	-	65cm	.038" (0.97mm)	25cm of Legato™ Hydrophilic Coating
410038VER-H	4	-	100cm	.038" (0.97mm)	40cm of Legato™ Hydrophilic Coating

NEW



CATHETERS

