



# Coaxial Introducers


8 | Accessories

9 | Mini Access Kits

May not be available for all OEM applications.

Single use. Sold packaged sterile. For non-sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS



Needles  
Page 11

## Accessories: Coaxial Introducer/Dilator Pair

Coaxial introducer/dilator pair - individually pouched; 5 units per box.

Catalog Number	Fench Size	Length	Description
<b>MAK4F-10</b>	4	10cm	
<b>MAK5F-10</b>	5	10cm	
<b>MAK5F-15</b>	5	15cm	
<b>S-MAK4F-10</b>	4	10cm	With stiffened dilator
<b>S-MAK5F-10</b>	5	10cm	With stiffened dilator
<b>S-MAK5F-15</b>	5	15cm	With stiffened dilator

## Accessories: Guide Wires

The wires are contained in a hoop & individually pouched; 10 units per box.

Catalog Number	Size	Description
<b>MAK001</b>	0.018" (0.46mm) x 40cm	Stainless steel guide wire with stainless steel tip
<b>MAK001N</b>	0.018" (0.46mm) x 40cm	Nitinol guide wire with platinum tip
<b>MAK001N-60</b>	0.018" (0.46mm) x 60cm	Nitinol guide wire with platinum tip
<b>MAK001N-80</b>	0.018" (0.46mm) x 80cm	Nitinol guide wire with platinum tip
<b>MAK005</b>	0.018" (0.46mm) x 40cm	Stainless steel guide wire with platinum tip
<b>MAK005-60</b>	0.018" (0.46mm) x 60cm	Stainless steel guide wire with platinum tip

## Accessories: Angiographic Needles

Packaged non-sterile, auto-bagged, 100 per case.

Catalog Number	Gauge	Length	Description
<b>AD21T21WP</b>	21G	2.5cm	
<b>AD21T41WP</b>	21G	4cm	
<b>AD21T71WP</b>	21G	7cm	
<b>AD21T91WP</b>	21G	9cm	
<b>AD21T41W-EP</b>	21G	4cm	Echo enhanced
<b>AD21T71W-EP</b>	21G	7cm	Echo enhanced
<b>AD21T91W-EP</b>	21G	9cm	Echo enhanced

# Merit MAK™ Mini Access Kits

The Merit MAK consists of a 4F or 5F coaxial introducer/dilator pair, 21 gauge introducer needle, and 0.018" (0.46mm) guide wire. 21G needles are available echo enhanced (4cm, 7cm & 9cm) and non-echo enhanced (7cm & 9cm). Needle and coaxial pair are pouched individually and contained in an outer kit pouch with the guide wire (guide wire is in a hoop). Both 4F and 5F products accept an 0.035" (0.89mm) or 0.038" (0.97mm) guide wire.

## STANDARD 7CM NEEDLE KITS

Catalog Number	French		Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
	Size	Length				
<b>MAK401</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
<b>MAK501</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
<b>MAK401N</b>	4	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
<b>MAK501N</b>	5	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
<b>MAK501N-15</b>	5	15cm	0.018" (0.46mm) x 60cm	NT/PT	21G x 7cm	10
<b>MAK405</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
<b>MAK505</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10


## ECHO ENHANCED 7CM NEEDLE KITS

Catalog Number	French		Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
	Size	Length				
<b>MAK401-E</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
<b>MAK501-E</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
<b>MAK401N-E</b>	4	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
<b>MAK501N-E</b>	5	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
<b>MAK405-E</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
<b>MAK505-E</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10

## ECHO ENHANCED 4CM NEEDLE KITS

Catalog Number	French		Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
	Size	Length				
<b>MAK401-PED</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	10
<b>MAK501-PED</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	10
<b>MAK401N-PED</b>	4	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	10
<b>MAK501N-PED</b>	5	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	10
<b>MAK405-PED</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 4cm	10
<b>MAK505-PED</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 4cm	10

SS = Stainless Steel    PT = Platinum Tip    NT = Nitinol



# Merit S-MAK™ Stiffened Dilator Mini Access Kits

The Merit S-MAK consists of a 4F or 5F coaxial introducer/stiffened dilator pair, 21 gauge echo enhanced introducer needle, and 0.018" (0.46mm) guide wire. Needle and coaxial pair are pouched individually and contained in an outer kit pouch with the guide wire (guide wire is in a hoop). Both 4F and 5F products accept an 0.035" (0.89mm) or 0.038" (0.97mm) guide wire.

## STIFFENED DILATOR KITS - Echo enhanced 7cm needle

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
<b>S-MAK401</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
<b>S-MAK501</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
<b>S-MAK401N</b>	4	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
<b>S-MAK501N</b>	5	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
<b>S-MAK501N-15</b>	5	15cm	0.018" (0.46mm) x 60cm	NT/PT	21G x 7cm	10
<b>S-MAK405</b>	4	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
<b>S-MAK505</b>	5	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10

SS = Stainless Steel    PT = Platinum Tip    NT = Nitinol

