

Guide Wire Accessories

GUIDE WIRES



Metal Guide Wire Insertion Tool
 Hub Material: Acrylic
 Inner Lumen: .022"
 Shaft Material: Stainless Steel

501994001



Plastic Guide Wire Insertion Tool
 Hub Material: High Density Polyethylene
 Inner Lumen: .028"
 Shaft Material: High Density Polyethylene

MAP600B



2,000cc Guide Wire Bowl

M-P2000



Ringmaster™ Guide Wire Basin

RM01



J-Straightener

101415003



J-Straightener

101415004



Natural Color Retention Clip

100730002