

Torque Devices (Pin Vise)

Cap Material: Polycarbonate Body Material: Polypropylene



Salmon. Fits Guide Wires
.009" to .018"

510061001

Green. Fits Guide Wires
.009" to .018"

510061002

Royal Blue. Fits Guide Wires
.009" to .018"

510061003

White. Fits Guide Wires
.009" to .018"

510061004

Salmon. Fits Guide Wires
.013" to .024"

502013001

Salmon. Fits Guide Wires
.021" to .040"

501985001

GUIDE WIRES

Sea Dragon™ Torque Device

H2O Torq™ Torque Device



Fits Guide Wires .018" (0.46mm) to
.038" (0.97mm) H2O

501695001



Orange. Fits Guide Wires
From .025" (0.64mm) to
.040" (1.02mm)

701655001



Yellow. Fits Guide Wires
From .010" (0.25mm) to
.020" (0.51mm)

701656001