

Hemostasis Valves & Y-Adapters

- 120 FLO 30™ Tuohy Borst
- 120 Access-9™
- 122 AccessPLUS™
- 122 DoublePlay™
- 122 Passage®
- 123 Honor®
- 123 MBA™
- 124 Hemostasis Valve Adapter

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.


SEE RELATED PRODUCTS



Guide Wires
Page 29



Catheters
Page 37



Needles
Page 11

FLO 30™ Tuohy Borst

Material: Polycarbonate Inner Lumen: 9F (.120") Luer: Fixed



Fixed Male Luer
French Size up to 9F

FLO30P

Access-9™

Material: Polycarbonate Inner Lumen: 9F (.120")
Valve Material: Silicone Luer: Rotating



Large Bore
Hemostasis Valve

500366002



Large Bore Hemostasis Valve
with 6" PVC Tubing extension

535001001

DEVICES



Large Bore Hemostasis Valve with
12" PVC Tubing Extension

535003012



Large Bore Hemostasis Valve with
24" Polyurethane Tubing Extension

535012024



Large Bore Hemostasis Valve with
12" Polyurethane Tubing Extension

535012012

AccessPLUS™

Material: Polycarbonate Inner Lumen: 9F (.120")
Valve Material: Silicone Luer: Rotating



Large Bore Hemostasis Valve
500628002



Large Bore Hemostasis Valve with
12" PVC Tubing Extension
563010001

DoublePlay™

Material: Polycarbonate
Valve Material: Silicone
Luer: Rotating



Large Bore Double Y Hemostasis Valve
Inner Lumen: 9F (.120")
501075001

Passage®

Material: Polycarbonate
Valve Material: Silicone



Passage® Hemostasis
Valve with Rotating Luer
500066



Hemostasis Valve with Slip Fit
(.179" to .209" taper)
Inner Lumen: 7.5F (.097")
K12-00457P

Honor®

Material: Polycarbonate Inner Lumen: 7F (.094")
Valve Material: Silicone Luer: Rotating



Hemostasis Valve
580195001



Hemostasis Valve with
8" Tubing Extension
MAP350

MBA™

Material: Polycarbonate Inner Lumen: 9F (.118")
Valve Material: Silicone Luer: Rotating



Large Bore Dual Seal
Hemostasis Valve
500981002



Hemostasis Valve with 10"
PVC Tubing Extension
501086001

Hemostasis Valve Adapter

Material: Polycarbonate Valve Material: Silicone
Cap: ABS Handle: Acetal

Adapter with Stopcock Handle
Male with Swivel Luer
French Size: 8F or Smaller

