


# Luer Activated Valves

- Proprietary Silicone Valve provides leak proof, swabable seal for needless luer access
- ISO 10993 Compliance: All materials are DEHP-free and latex free
- Polycarbonate and copolyester materials for easy bonding

Merit's Luer Activated Valves are designed to aspirate or inject fluids on demand. The proprietary valve facilitates multiple uses and, upon disconnection, does not require a dead-end cap to prevent leaking.

SEE RELATED PRODUCTS



Syringes  
Page 225



Stopcocks  
Page 199

# Swabable, Needleless Valves

Body: Polycarbonate

Stem: Silicone



Y-Port 4.0mm, PC  
**702238001**



Male Luer Lock, PC  
**701000001**



Male Luer Lock, Red, PC  
**702242001**



Male Luer Lock, Copolyester  
**702237001**



Y-Port 4.0mm, Copolyester  
**702239001**



T-Port 6.8mm ID Valve Adaptor  
**702236001**



Tube End 2.0mm, PC  
**702243001**

## Swabable, Needleless Valves



Tube End 2.5mm, PC

**702243002**



Tube End 2.8mm, PC

**702243003**



Tube End 3.7mm, PC

**702243004**



Tube End 4.0mm, PC

**701814001**



Bag Valve 6.8mm OD, PC

**702240001**



Vial Adapter 13mm, PC

**702244001**



Vial Adapter 20mm, PC

**702244002**



Breakoff Valve 6.0mm, PC

**702245001**



Breakoff Valve 6.6mm, PC

**702245002**



Male Slip Luer, PC

**702241001**