


Luers & Connectors

- Male and female
- Slip luer
- Tubing pockets

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS



Cap & Covers
Page 201



Tubing
Page 191

Female Slip Fit Luers

Material: Polycarbonate



White Female Luer Lock
Slip Fit .120 to .150

100916001



Clear Female Luer Lock
Slip Fit .120 to .150 Taper

102112001



Clear Female Luer Lock
Slip Fit .154 to .192 Taper

102102001



Gamma Color Female Luer Lock
Slip Fit .154 to .192 Taper

102202001



Gamma Color Female Luer Lock
Slip Fit .120 to .150 Taper

102202002



Clear Female Luer Lock
Slip Fit .179 to .209

100719

Female Tubing Pocket Fit Luers

Material: Polycarbonate



Clear Female Luer Lock
Tubing Pocket .095 to .101

102103002



Royal Blue Female Luer Lock
Tubing Pocket .095 to .101

100972002



Clear Female Luer Lock
Tubing Pocket .118 to .124

101023001



Yellow Female Luer Lock
Tubing Pocket .118 to .124

101023002



Clear Female Luer Lock
Tubing Pocket .124 to .132

102103001



Clear Female Luer Lock
Tubing Pocket .140 to .145

102103004



Clear Female Luer Lock
Tubing Pocket .158 to .164

102103003



Clear Female Luer Lock
Tubing Pocket .182 to .192

100976001



Gamma Color Female Luer Lock
Tubing Pocket .095 to .101

102203002



Gamma Color Female Luer Lock
Tubing Pocket .118 to .124

101123001



Gamma Color Female Luer Lock
Tubing Pocket .124 to .132

102203001



Gamma Color Female Luer Lock
Tubing Pocket .140 to .145

102203003



Clear Female Luer Lock
Tubing Pocket .124 to .132

101071002

Male Slip Fit Luers

Material: Polycarbonate



Blue Male Luer Lock
Slip Fit .120 to .150 Taper
100914002



Gamma Color Male Luer Lock
Slip Fit .120 to .150 Taper
102200002



Gamma Color Male Luer Lock
with Slip Fit .154 to .192
102200001



White Male Luer Lock
Slip Fit .120 to .150 Taper
100914001



Clear Male Luer Lock
Slip Fit .120 to .150 Taper
102110001



Clear Male Luer Lock
Slip Fit .154 to .192
102100001



Gamma Color, Male Luer
Lock Slip Fit, .148 to .166
102056002



Gamma Color Male Luer Lock
Slip Fit .178 to .208
102039002



Clear Male Luer Lock,
Slip Fit, .148 to .166
102056001



Clear Male Luer Lock
Slip Fit .178 to .208
102039

Male Tubing Pocket Fit Luers

Material: Polycarbonate



Clear Male Luer Lock
Tubing Pocket .095 to .101

102101002



Clear Male Luer Lock
Tubing Pocket .108 to .114

102101005



Clear Male Luer Lock
Tubing Pocket .118 to .124

101024001



Clear Male Luer Lock
Tubing Pocket .124 to .132

102101001



Clear Male Luer Lock
Tubing Pocket .140 to .145

102101004



Clear Male Luer Lock
Tubing Pocket .158 to .164

102101003



Clear Male Luer Lock
Tubing Pocket .182 to .192

100975001



Gamma Color Male Luer Lock
Tubing Pocket .118 to .124

101124001



Gamma Color Male Luer Lock
Tubing Pocket .124 to .132

102201001



Gamma Color Male Luer Lock
Tubing Pocket .095 to .101

102201002



Gamma Color Male Luer Lock
Tubing Pocket .158 to .164

102201003

Luers & Connectors Other



Clear, 2-Way Connector
Tubing Pockets
Connector 1: 0.160" Tubing Pocket
Connector 2: 0.160" Tubing Pocket
Connector 3: 0.160" Tubing Pocket
Material: Cyrolite

101002



Clear, 3-Way Connector
Tubing Pockets
Connector 1: 0.160" Tubing Pocket
Connector 2: 0.160" Tubing Pocket
Connector 3: 0.160" Tubing Pocket
Material: Cyrolite

101003



Y-Connector
Material: Polycarbonate

700155002



Clear Male/Female Luer
Lock Connector
Connector 1: Male Luer Lock
Connector 2: Female Luer Lock
Material: Polycarbonate

101047001



Clear Male/Male Luer
Lock Connector
Connector 1: Male Luer Lock
Connector 2: Male Luer Lock
Material: Polycarbonate

101025001



Clear Female/Female Luer
Lock Connector
Connector 1: Female Luer Lock
Connector 2: Female Luer Lock
Material: MABS

700357001



Clear, 2-Way Connector
Tubing Pockets
Connector 1: .0160" Tubing Pocket
Connector 2: 0.160" Tubing Pocket
Connector 3: 0.160" Tubing Pocket
Material: Polycarbonate

101004



Y Tubing Connector
Connector 1: .154 slip
Connector 2: .154 slip
Connector 3: .154 slip
Material: Clear ABS

101508001

T-Connectors

Material: Polycarbonate



Clear Two Female Luer Lock,
One Male Slip

100966001



Gamma Color Two Female
Luer Lock, One Male Slip

100966002



Clear One Female Luer Lock,
Two Male Slip

100969001



Clear All Male Slip

100967001



Clear Two Female Luer Lock,
One Male Swivel

501314001



Clear All Female Luer Lock

100965001



Clear with 0.139" Tubing Pockets

100965002