

Marker Bands

CUSTOM

BUILT TO ORDER

In 2007, Merit acquired a technology from Lightek Corporation to create flexible radiopaque plastic laminate composition marker bands.

This technology enables Merit to place polyurethane marker bands, impregnated with Tungsten, Platinum, or other precious metal powder, on polyurethane tubing.

The marker bands are tough, flexible, and visible under fluoroscopy.

- Catheter shaft OD's down to 0.020"
- Wall thickness down to 0.003"
- Ideal for urethane family catheters



Platinum
Marker
Band



Used in Merit's own
Drainage Catheter Line