

Catheters

Microcatheters


May not be available for all OEM applications.

Single use. Sold packaged sterile. For non-sterile and other packaging options, please contact a Merit OEM representative.


SEE RELATED PRODUCTS



Hemostasis
Valves



Guide Wires



Needles

Microcatheters


A multipurpose microcatheter designed for use in peripheral and coronary vasculature.
Available in two sizes to accommodate most coils, embolic particles and sphericals.

Flexible distal 200mm designed for atraumatic vessel entry.

Hydrophilic coating (distal 80cm) for smooth traction through tortuous anatomy.

New swan neck design high performance of preshaped distal end ideal for access of hepatic artery.

Radiopaque Marker for easy identification of distal tip of catheter. Located 1.3mm proximal to the distal tip.



Why Use Merit Maestro® Microcatheters?

- All Merit Maestro Microcatheters are rated to 600psi
- (2) 3ml Medallion® syringes included for ease in mixing embolics/sphericals
- Embolics up to 700um (microns)
- Coils .018"

Merit Maestro Microcatheters 2.8F Tapered to 2.4F

Catalog Number	Length		Tip	Minimum Inner Diameter	Guide Wire	Guide Catheter Indicated for Guide Wire Size
	Proximal	Distal				
28MC24110ST	2.8F	2.4F	110	Straight	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC24130ST	2.8F	2.4F	130	Straight	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC24150ST	2.8F	2.4F	150	Straight	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC2411045	2.8F	2.4F	110	45 degree	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC2413045	2.8F	2.4F	130	45 degree	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC2415045	2.8F	2.4F	150	45 degree	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC24110SN	2.8F	2.4F	110	Swan Neck	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC24130SN	2.8F	2.4F	130	Swan Neck	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC24150SN	2.8F	2.4F	150	Swan Neck	.020" (0.53mm)	.018" (0.46mm) .038" (0.97mm)
28MC28150SN	2.8F	2.8F	150	Swan Neck	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)

Merit Maestro Microcatheters 2.8 No Taper

Catalog Number	Length		Tip	Minimum Inner Diameter	Guide Wire	Guide Catheter Indicated for Guide Wire Size
	Proximal	Distal				
28MC28110ST	2.8F	2.8F	110	Straight	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC28130ST	2.8F	2.8F	130	Straight	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC28150ST	2.8F	2.8F	150	Straight	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC2811045	2.8F	2.8F	110	45 degree	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC2813045	2.8F	2.8F	130	45 degree	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC2815045	2.8F	2.8F	150	45 degree	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC28110SN	2.8F	2.8F	110	Swan Neck	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC28130SN	2.8F	2.8F	130	Swan Neck	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)
28MC28150SN	2.8F	2.8F	150	Swan Neck	.024" (0.64mm)	.021" (0.53mm) .038" (0.97mm)