

# Safety & Waste Collection

Disposal Basins

Drainage Bags

Fluid & Waste Disposal Systems

Sharps Holders

Pen & Labels


See also SecureLoc™ Needles

See also Futura® Safety Scalpels


May not be available for all OEM applications.

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.


SEE RELATED PRODUCTS



Syringes



Scalpels



Needles

## Disposal Basins



BackStop® Disposal Basin, 500ml fluid capacity and absorbent material. Includes particulate-free die cut sponge and seal-tight lid.

**MBS100P**



BackStop+™ Disposal Basin with temporary sharps holder, 500ml fluid capacity and absorbent material. Includes particulate-free die cut sponge and seal-tight lid.

**MBSP100P**



DugOut® Disposal Basin 1000ml fluid capacity and absorbent material. Includes particulate-free die cut sponge, secondary compartment for storage of ancillary products or saline and seal-tight lid.

**DUG100P**



DugOut Disposal Basin with 1000ml fluid capacity and absorbent material. Includes particulate-free die cut sponge, secondary compartment for storage of ancillary products or saline and seal-tight lid with ShortStop™.

**DUG1SSP**

## Disposal Basins



MiniStop+™ Disposal Basin with temporary sharps holder, 100ml fluid capacity and absorbent material. Includes particulate-free die cut sponge and seal-tight lid.

**MMSP100P**



MiniStop® Disposal Basin, 100ml fluid capacity and absorbent material. Includes particulate-free die cut sponge and seal-tight lid.

**MMS100P**

## Drainage Bags



Merit Drainage Depot™ 600ml Nephrostomy Bag, 24" (61cm). Kink-resistant Adjustable Tubing with Fixed Male Luer, Twist Drain Valve, Velcro® Brand Fastener and Additional Male Luer.

**MDD600P**



Merit Drainage Depot™ 600ml Nephrostomy Bag with Stepped Connector, Twist Drain Valve and Velcro® Brand Fastener.

**MDD610P**

## Drainage Bags



600ml bag, soft cloth backing, stepped connector, twist drain valve and VELCRO® brand fastener.

**FZ610P**



600ml bag, soft cloth backing, 24" (61cm) kink-resistant adjustable tubing with fixed male luer, twist drain valve and VELCRO® brand fastener. Additional male luer.

**FZ624P**



1000ml bag, soft cloth backing, 24" (61cm) kink-resistant adjustable tubing with fixed male luer, twist drain valve and VELCRO® brand fastener. Additional male luer.

**FZ1000P**

## Fluid & Waste Disposal Systems



Merit Disposal Depot, 1400ml Bag,  
72" (183cm) Tubing, Fixed Male Luer

**MDD100P**



Merit Disposal Depot, 1400ml Bag, 72" (183cm)  
Tubing, One-Way Check Valve with Female Luer

**MDD110P**



Merit Disposal Depot, 1400ml bag, Check Relief  
Valve to Vented Spike, Female Control Port  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**MDD200HP**



Merit Disposal Depot, 1400ml Bag, Check Relief  
Valve to Vented Spike, Female Control Port with Reflux Valve  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**MDD220HP**

## Fluid & Waste Disposal Systems



Merit Disposal Depot, 1400ml Bag,  
Check Relief Valve to Vented Spike,  
Rotating Male Connector Control Port  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**MDD300HP**



Merit Disposal Depot, 1400ml bag,  
Check Relief Valve to Vented Drip Chamber,  
Rotating Male Connector Control Port  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**MDD400HP**



Merit Disposal Depot, 1400ml Bag, Check Relief  
Valve to Vented Spike, 3-Way Stopcock with  
Vented Spike, Female Control Port  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**MDD500HP**



Merit Disposal Depot™, 1400ml Bag, Check Relief  
Valve to Vented Spike, 3-Way Stopcock with  
Vented Spike, Female Control Port with Reflux Valve  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**MDD520HP**

## Fluid & Waste Disposal Systems



Merit Disposal Depot, 1400ml Bag, Check Relief Valve to Vented Spike, 3-Way Stopcock with In-Line Contrast Management System, Female Control Port with Reflux Valve  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**K10-00577P**



Merit Disposal Depot, 1400ml Bag, Check Relief Valve to Vented Spike, 3-Way Stopcock with Miser Contrast Management System, Female Control Port with Reflux Valve  
Valve to Bag: 48" (122cm)  
Valve to Spike/Drip: 72" (183cm)

**K10-00576P**



Captiva® Blood  
Containment Device

**CAN100P**

## Sharps Holders



ShortStop® Temporary Sharps Holder

**SS10P**



ShortStop Advantage Temporary Sharps Holder with Side Access Point

**SS20P**





Shortstop Grandstand Temporary Sharps Holder with Upright Access

**SS30P**

## Pen & Labels

Packaging: Sterile Single Pouch - 25/Box and 100/Case



Catalog Number

Description

- PAL100** Contains pen and label set of 35
- PAL200** Contains pen and label set of 24
- PAL225** Contains pen, label set of 24 and guide wire bowl label set of 5
- PAL250** Contains pen, label set of 24 and large bowl label set of 7
- PAL300** Contains pen and label set of 12
- PAL325** Contains pen, label set of 12 and guide wire bowl label set of 5
- PAL350** Contains pen, label set of 12 and large bowl label set of 7
- PAL400** Contains pen and blank label set of 12
- PAL425** Contains pen, blank label set of 12 and guide wire bowl label set of 5
- PAL450** Contains pen, blank label set of 12 and large bowl label set of 7
- PAL500** Contains pen and guide wire bowl label set of 5
- PAL600** Contains pen and large bowl label set of 7
- PAL700** Contains pen, label set of 35 and large bowl label set of 7
- K13-00476** Pen only