

Scalpels

Futura® retractable safety scalpels are designed to protect users from accidental cuts while passing and disposing.

- Blade sizes are color-coded for easy identification
- Small, compact and ideal for kits, trays and packs
- Made of latex free materials

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS




Sharps Holders
Page 168



Needles
Page 11

Futura® Safety Scalpels

Material: Polystyrene Blade Material: 440 A Steel



Blade Style #11, Blue
SMS211P



Blade Style #15, Yellow
SMS215P



Blade Style #10, Green
SMS210P