

# Marquis® Stopcocks

- 224 | 1-Way
- 225 | 3-Way
- 227 | 4-Way
- 228 | Customization

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.


## SEE RELATED PRODUCTS



Tubing  
Page 191



Custom Kits  
Page 141



Caps & Covers  
Page 201

# 1-Way Stopcocks

Large Bore: N Handle: Acetal Material: Polycarbonate

1200 PSI



Fixed Male Luer  
**U1FLP**



Swivel Nut  
**U1SNP**



Right Rotating Male Luer  
**U1RP**



Slip Luer  
**U1SLP**

1050 PSI



Fixed Male Luer  
**H1FLP**



Swivel Nut  
**H1SNP**



Right Rotating Male Luer  
**H1RP**

500 PSI



Fixed Male Luer  
**M1FLP**



Swivel Nut  
**M1SNP**



Slip Luer  
**M1SLP**

200 PSI



Fixed Male Luer  
**S1FLP**



Swivel Nut  
**S1SNP**



Slip Luer  
**S1SLP**

# 3-Way Stopcocks

Handle: Acetal Material: Polycarbonate

## 1200 PSI



Slip Luer  
**U3SLP**



Fixed Male Luer  
**U3FLP**



Right Rotating Male Luer  
**U3RRP**



Swivel Nut  
**U3SNP**

## 1050 PSI



Slip Luer  
**H3SLP**



Fixed Male Luer  
**H3FLP**



Left Rotating Male Luer  
**H3RLP**



Right Rotating Male Luer  
**H3RRP**



Swivel Nut  
**H3SNP**

## 500 PSI



Slip Luer  
**M3SLP**



Fixed Male Luer  
**M3FLP**



Left Rotating Male Luer  
**M3RLP**



Right Rotating Male Luer  
**M3RRP**



Swivel Nut  
**M3SNP**

# 3-Way Stopcocks

Handle: Acetal Material: Polycarbonate



Swivel Nut  
**S3SNP**



Slip Luer  
**S3SLP**



All Female Luer  
**S3FP**



Fixed Male Luer  
**S3FLP**

200 PSI



Barbed Connector,  
Red Handle  
**501641004**



Barbed Connector,  
Green Handle  
**501641006**



Barbed Connector,  
Blue Handle  
**501641008**



Swivel Nut, Large Bore  
**L3SNP**



With Tubing Pocket .158  
**501868001**



With Tubing Pocket .110  
**501855001**

50 PSI

# 4-Way Stopcocks

Handle: Acetal Material: Polycarbonate

1200 PSI



Right Rotating Male Luer  
**U4RRP**

Swivel Nut  
**U4SNP**

Fixed Male Luer  
**U4FLP**

Slip Luer  
**U4SLP**

1050 PSI



Fixed Male Luer  
**H4FLP**

Right Rotating Male Luer  
**H4RRP**

Swivel Nut  
**H4SNP**

500 PSI



Right Rotating Male Luer  
**M4RRP**

Swivel Nut  
**M4SNP**

Fixed Male Luer  
**M4FLP**

Right Rotating Male Luer w/ On Handle  
**MN4RRP**

# 4-Way Stopcocks

Handle: Acetal Material: Polycarbonate



50 PSI, Swivel Nut,  
**Large Bore**  
**L4SNP**



200 PSI,  
Swivel Nut  
**S4SNP**



200 PSI,  
Fixed Male Luer  
**S4FLP**



200 PSI,  
Slip Luer  
**S4SLP**

50 to 200 PSI

## Customization Examples



Contact a Merit Representative about customization options including handle color, port configuration, tubing pockets and assembly.