


# Syringes

- 230 X-Change™
- 231 Coronary Control
- 237 Medallion®
- 247 Medallion® COP
- 248 VacLok®

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS



Caps & Covers  
Page 201



Inflation Devices  
Page 125



Stopcocks  
Page 223

# Merit X-Change™ Syringe

Plunger Material: ABS

Material: Polycarbonate

Tip Material: Silicone

**NEW**



20ml open-ended barrel syringe  
with removable fixed male luer tip

**EX20-MLP**

20ml open-ended barrel syringe  
with removable slip luer tip

**EX20-SLP**

Dead-End Cap

**102218001**

## CUSTOMIZING

From custom pad printing on the barrel, custom colored plungers, or custom barrel tip configurations, Merit Medical can work with your product design team to customize the X-Change to fit your needs.

# Coronary Control Syringes

With Ring Plunger and Ring Grip, Rotating Adapter

Plunger Material: ABS Material: Polycarbonate

Tip Material: Silicone



8ml, Inject8™, Solid Plunger

**CCS880P**



10ml, Inject10n™, Solid Plunger

**CCX010P**



10ml, Solid Plunger

**CCS800P**



12ml, Solid Plunger

**CCS100P**



12ml, Solid Plunger

**CCS600P**



20ml, Solid Plunger

**CCS320P**

# Coronary Control Syringes

With Ring Plunger and Ring Grip, Fixed Male Connector

Plunger Material: ABS Material: Polycarbonate  
Tip Material: Silicone



8ml, Inject8, Solid Plunger  
**CCS881P**



10ml, Inject10n, Solid Plunger  
**CCX011P**



10ml, Solid Plunger  
**CCS801P**



12ml, Solid Plunger  
**CCS601P**



20ml, Solid Plunger  
**CCS321P**



10ml, Open Bore  
**K06-00189P**



20ml, Open Bore  
**K06-00281P**

## Open Bore Control Syringes

With Palm Pad Plunger and Ring Grip, Open Bore

## Coronary Control Smart Tip™ Syringes

Rotating Adapter

Plunger Material: ABS Material: Polycarbonate

Tip Material: Blue Silicone



10ml, Smart Tip, Solid Plunger

**ST800P**



12ml, Smart Tip, Solid Plunger

**ST600P**



20ml, Smart Tip, Solid Plunger

**ST320P**

## Coronary Control Smart Tip Syringes

Fixed Male Connector

Plunger Material: ABS Material: Polycarbonate

Tip Material: Blue Silicone



10ml, Smart Tip, Solid Plunger

**ST801P**



12ml, Smart Tip, Solid Plunger

**ST601P**



20ml, Smart Tip, Solid Plunger

**ST321P**

# Coronary Control Syringes

## Locking Ring Plunger

Plunger Material: ABS Material: Polycarbonate  
Tip Material: Silicone



6ml, Locking Ring Plunger  
Fixed Male Adapter

**CCS461P**



6ml, Locking Ring Plunger  
Rotating Male Adapter

**CCS460P**



12ml, Locking Ring Plunger  
Fixed Male Adapter

**CCS101P**



12ml, Locking Ring Plunger  
Rotating Male Adapter

**CCS200P**



12ml, Locking Ring Plunger  
Fixed Male Adapter

**CCS201P**

# Coronary Control Syringes

With Palm Pad Plunger and Ring Grip, Rotating Adapter

Plunger Material: ABS Material: Polycarbonate

Tip Material: Silicone



8ml, Inject8,  
Solid Plunger

**CCSB880P**



10ml, Inject10n,  
Solid Plunger

**CCXB010P**



10ml, Solid Plunger

**CCSB900P**



12ml, Solid Plunger

**CCSB700P**

# Coronary Control Syringes

With Palm Pad Plunger and Ring Grip, Fixed Male Connector

Plunger Material: ABS Material: Polycarbonate

Tip Material: Silicone



8ml, Inject8,  
Solid Plunger

**CCSB881P**



10ml, Inject10n,  
Solid Plunger

**CCXB011P**



10ml, Solid Plunger

**CCSB901P**



12ml, Solid Plunger

**CCSB701P**

## Coronary Control Syringes

With Palm Pad Plunger and Wing Grip, Rotating Adapter

Plunger Material: ABS Material: Polycarbonate  
Tip Material: Silicone



8ml, Inject8,  
Solid Plunger  
**CCSW880P**



10ml, Solid Plunger  
**CCSB910P**



12ml, Solid Plunger  
**CCSB710P**

## Coronary Control Syringes

With Palm Pad Plunger and Wing Grip, Fixed Male Connector

Plunger Material: ABS Material: Polycarbonate  
Tip Material: Silicone



8ml, Inject8,  
Solid Plunger  
**CCSW881P**




10ml, Solid Plunger  
**CCSB911P**



12ml, Solid Plunger  
**CCSB711P**

# Medallion® Syringes



Note: Slip Luer connectors are available upon special request.

## 1ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



White

**MSS011P**

Light Blue

**MSS011-LBP**

Yellow

**MSS011-YP**

Red

**MSS011-RP**

Purple

**MSS011-PRP**

Light Green

**MSS011-LGP**

Dark Green

**MSS011-DGP**

## 3ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



White

**MSS031P**

Light Blue

**MSS031-LBP**

Yellow

**MSS031-YP**

Red

**MSS031-RP**

Purple

**MSS031-PRP**

Light Green

**MSS031-LGP**

Dark Green

**MSS031-DGP**

## 6ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



White

**MSS061P**

Light Blue

**MSS061-LBP**

Yellow

**MSS061-YP**

Red

**MSS061-RP**

Purple

**MSS061-PRP**

Light Green

**MSS061-LGP**

Dark Green

**MSS061-DGP**

## 10ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



White

**MSS111P**

Light Blue

**MSS111-LBP**

Yellow

**MSS111-YP**

Red

**MSS111-RP**

Purple

**MSS111-PRP**

Light Green

**MSS111-LGP**

Dark Green

**MSS111-DGP**

## 20ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



White

**MSS121P**

Light Blue

**MSS121-LBP**

Yellow

**MSS121-YP**

Red

**MSS121-RP**

Purple

**MSS121-PRP**

Light Green

**MSS121-LGP**

Dark Green

**MSS121-DGP**

## 30ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



White

**MSS131P**

Light Blue

**MSS131-LBP**

Yellow

**MSS131-YP**

Red

**MSS131-RP**

Purple

**MSS131-PRP**

Light Green

**MSS131-LGP**

Dark Green

**MSS131-DGP**

## 60ml Fixed Male Medallion Syringes

Barrel Material: Polycarbonate  
Tip Material: Silicone  
Plunger Material: ABS

## Customized Medallion Syringes

Barrel Material: Polycarbonate  
Tip Material: Silicone  
Plunger Material: ABS



White

**MSS160P**

Red

**MSS160-RP**

Light Blue

**MSS160-LBP**

10ml White  
without Logo

**K01-05586P**

6ml White  
without Logo

**K01-07201P**

6ml Yellow  
without Logo

**K01-06124P**

20ml Green with  
Gamma Barrel

**K01-06182P**

# 10ml Sword Handle Fixed Male Medallion Syringes

Barrel Material: Polycarbonate

Tip Material: Silicone

Plunger Material: ABS



White

**MSSW11P**

Light Blue

**MSSW11-LBP**

Yellow

**MSSW11-YP**

Red

**MSSW11-RP**

Purple

**MSSW11-PRP**

Light Green

**MSSW11-LGP**

Dark Green

**MSSW11-DGP**

## .25ml Sword Handle Fixed Male Medallion Syringe

Barrel Material: Polycarbonate

Tip Material: Silicone

Plunger Material: ABS

## 1ml Polypropylene Medallion Syringe

Barrel Material: Polypropylene

Tip Material: Silicone

Plunger Material: Acrylic Copolymer



White

**K01-05984P**



**K01-05934P**

## 20ml Sword Handle Fixed Male Medallion Syringes

Barrel Material: Polycarbonate  
Tip Material: Silicone  
Plunger Material: ABS



White

**MSSW21P**

Light Blue

**MSSW21-LBP**

Yellow

**MSSW21-YP**

Red

**MSSW21-RP**

Purple

**MSSW21-PRP**

Light Green

**MSSW21-LGP**

Dark Green

**MSSW21-DGP**

## Medallion Slip Luer Syringes

Barrel Material: Polycarbonate  
Tip Material: Silicone  
Plunger Material: ABS



3ml, Dark Green

**K01-02040P**

6ml, Yellow

**K01-02044P**

10ml, Dark Green

**K01-02025P**

10ml, Yellow

**K01-60131P**

1ml, Light Blue

**K01-07287P**

1ml, White

**K01-05235P**

## Toomey Tip Medallion Syringe



60ml, Toomey tip  
**MSS160TP**

## RXP® Flush Syringe



4ml, Rapid Exchange Prep Syringe  
**RXPORP**

## NITRO Labeled, 10ml Fixed Male Luer

Barrel Material: Polycarbonate  
Tip Material: Silicone Plunger Material: ABS



White

**K01-60077P**

Light Blue

**K01-60085P**

Yellow

**K01-60101P**

Red

**K01-60109P**

Light Green

**K01-60093P**

## LIDOCAINE Labeled, 10ml

Fixed Male Luer

Barrel Material: Polycarbonate  
Tip Material: Silicone Plunger Material: ABS



Light Blue

**K01-60081P**

Yellow

**K01-60097P**

Red

**K01-60105P**

Light Green

**K01-60089P**

Purple

**K01-04086P**

## HEPARIN Labeled, 10ml

Fixed Male Luer

Barrel Material: Polycarbonate  
Tip Material: Silicone Plunger Material: ABS



Light Blue

**K01-60082P**

Yellow

**K01-60098P**

Red

**K01-60106P**

Light Green

**K01-60090P**

# CONTRAST Labeled, 10ml

Fixed Male Luer

Barrel Material: Polycarbonate  
Tip Material: Silicone Plunger Material: ABS



White	Light Blue	Yellow	Red	Light Green	Dark Green
<b>K01-60078P</b>	<b>K01-60086P</b>	<b>K01-60102P</b>	<b>K01-60110P</b>	<b>K01-60094P</b>	<b>K01-00786P</b>

# SALINE Labeled, 10ml and 20ml

Medication Printed, Fixed Male Luer

Barrel Material: Polycarbonate  
Tip Material: Silicone  
Plunger Material: ABS

# VERSED Labeled, 3ml

Medication Printed, Fixed Male Luer

Barrel Material: Polycarbonate  
Tip Material: Silicone  
Plunger Material: ABS



White, 10ml	Light Blue, 10ml	Yellow, 10ml	Light Blue, 20ml, sword handle	Light Blue
<b>K01-60075P</b>	<b>K01-60083P</b>	<b>K01-60099P</b>	<b>K01-6099P</b>	<b>K01-05185P</b>

# XYLOCAINE® Labeled, 10ml

Medication Printed, Fixed Male Luer

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



Light Blue  
**K01-00917P**



Red  
**K01-00953P**

Xylocaine® is a registered trademark of APP Pharmaceuticals, LLC.

# FLUSH Labeled, 10ml

Medication Printed, Fixed Male Luer

Barrel Material: Polycarbonate Tip Material: Silicone  
Plunger Material: ABS



Light Blue  
**K01-05093P**



White  
**K01-05025P**

# Medallion® COP Syringes

## 20ml Fixed Male Medallion COP Syringes

Barrel Material: Cyclo Olefin Polymer (COP)

Tip Material: Silicone

Plunger Material: Polycarbonate or ABS



Red  
Polycarbonate  
Plunger

**COP121PC-RP**

Green  
Polycarbonate  
Plunger

**COP121PC-GP**

Purple  
Polycarbonate  
Plunger

**COP121PC-PP**

Yellow  
Polycarbonate  
Plunger

**COP121PC-YP**

Orange  
Polycarbonate  
Plunger

**COP121PC-OP**

Blue  
Polycarbonate  
Plunger

**COP121PC-BP**

Gray  
Polycarbonate  
Plunger

**COP121PC-GRP**



The Medallion® COP syringe is shipped in a protective tray.



Red  
ABS Plunger

**COP121AB-RP**

Green  
ABS Plunger

**COP121AB-GP**

Purple  
ABS Plunger

**COP121AB-PP**

Yellow  
ABS Plunger

**COP121AB-YP**

# VacLok® Negative Pressure Syringes

Plunger Material: ABS Material: Polycarbonate  
Tip Material: Silicone Luer: Fixed Male

## Four Position Locking Syringe



VacLok 10ml  
Negative Pressure  
White

**VAC110P**



VacLok 20ml  
Negative Pressure  
White

**VAC120P**



VacLok 20ml  
Negative Pressure  
Red

**VAC120-RP**



VacLok 20ml  
Negative Pressure  
Yellow

**VAC120-YP**



VacLok 20ml  
Negative Pressure  
Light Blue

**VAC120-LBP**

## Six Position Locking Syringe



VacLok 30ml  
Negative Pressure  
White

**VAC130P**



VacLok 30ml  
Negative Pressure  
Red

**VAC130-RP**



VacLok 30ml  
Negative Pressure  
Yellow

**VAC130-YP**



VacLok 30ml  
Negative Pressure  
Light Blue

**VAC130-LBP**



VacLok 60ml  
Negative Pressure  
White

**VAC160P**

# FEATURED PRODUCT

## VacLok® Locking Vacuum Syringe



### SPECIFICATIONS



Plunger Material:  
ABS

Barrel Material:  
Polycarbonate

Plunger Tip  
Material: Silicone  
Luer: Fixed Male

### WHAT IT DOES

Merit Medical's VacLok® is a single-use syringe designed to create a vacuum when withdrawal of a viscous material or evacuation of a device such as a balloon during a procedure is required.

### FEATURES & BENEFITS

Feature	Benefit
Patented Design	Locking plunger design eliminates need to manually maintain a vacuum
Sturdy Design	Polycarbonate barrel and ABS plunger provides a stable platform
Large Wings	Easy to grip
Multiple Locking Positions	Allows for multiple vacuum and volume settings

### CUSTOMIZING

From custom pad printing on the barrel to custom colored plungers, Merit Medical can work with your product design team to customize the VacLok to fit your needs. For sterile and other packaging options, please contact a Merit Medical OEM representative.

