

# Tubing

190 | High Pressure

192 | Pressure Monitoring

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS



Stopcocks  
Page 223



Check Relief  
Valves  
Page 137




Caps and  
Covers  
Page 201

# High Pressure Tubing

Rated at 1200 PSI (83 BAR).

Material: Flexible Nylon Braided Polyurethane

Male to Female



## Airless Rotating Adapter

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	72" Length (183cm)
<b>HPF100P</b>	<b>HPF200P</b>	<b>HPF300P</b>	<b>HPF480P</b>	<b>HPF720P</b>

## Fixed Male Connector


10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	72" Length (183cm)
<b>HPF101P</b>	<b>HPF201P</b>	<b>HPF301P</b>	<b>HPF481P</b>	<b>HPF721P</b>

# High Pressure Tubing

Rated at 1200, 900 and 500 PSI (83, 62 and 34 BAR).

Material: PVC, Clear

Male to Female



## Airless Rotating Adapter

1200 PSI (83 BAR)  
900 PSI (62 BAR)  
500 PSI (34 BAR)

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	59" Length (150cm)
<b>HPT100P</b>	<b>HPT200P</b>	<b>HPT300P</b>	<b>HPT480P</b>	<b>HPT590P</b>
<b>HP9100P</b>	<b>HP9200P</b>	<b>HP9300P</b>	<b>HP9480P</b>	
<b>HP5100P</b>	<b>HP5200P</b>	<b>HP5300P</b>	<b>HP5480P</b>	

## Fixed Male Connector

1200 PSI (83 BAR)  
900 PSI (62 BAR)  
500 PSI (34 BAR)

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	59" Length (150cm)
<b>HPT101P</b>	<b>HPT201P</b>	<b>HPT301P</b>	<b>HPT481P</b>	<b>HPT591P</b>
<b>HP9101P</b>	<b>HP9201P</b>	<b>HP9301P</b>	<b>HP9481P</b>	
<b>HP5101P</b>	<b>HP5201P</b>	<b>HP5301P</b>	<b>HP5481P</b>	







# Excite™ High Pressure Tubing

Rated at 1200 PSI (83 BAR).

Material: Coextruded Nylon/Polyurethane

Male to Female

## Airless Rotating Adapter

					
10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	39" Length (100cm)	48" Length (122cm)	59" Length (150cm)
<b>HPC100P</b>	<b>HPC200P</b>	<b>HPC300P</b>	<b>HPC390P</b>	<b>HPC480P</b>	<b>HPC590P</b>

## Fixed Male Connector








10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	39" Length (100cm)	48" Length (122cm)	59" Length (150cm)
<b>HPC101P</b>	<b>HPC201P</b>	<b>HPC301P</b>	<b>HPC391P</b>	<b>HPC481P</b>	<b>HPC591P</b>












# Pressure Monitoring Tubing

High Durometer, Non Expanding.

Material: PVC, Clear

							
<b>Male to Male</b>	6" Length (15 cm)	12" Length (30 cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	60" Length (152cm)	72" Length (183cm)
	<b>PM6006P</b>	<b>PM6012P</b>	<b>PM6024P</b>	<b>PM6036P</b>	<b>PM6048P</b>	<b>PM6060P</b>	<b>PM6072P</b>

						
<b>Male to Female</b>	6" Length (15 cm)	12" Length (30 cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	60" Length (152cm)
	<b>PM6106P</b>	<b>PM6112P</b>	<b>PM6124P</b>	<b>PM6136P</b>	<b>PM6148P</b>	<b>PM6160P</b>

			
<b>Male to Female</b>	72" Length (183cm)	84" Length (213 cm)	96" Length (244cm)
	<b>PM6172P</b>	<b>PM6184P</b>	<b>PM6196P</b>
















**NEW**

# Merit Freedom™ Pressure Monitoring Tubing

Material: PVC, Clear non-DEHP

Rated at 1200 PSI (83 BAR). High Durometer, Non Expanding.

						
Male to Male	6" Length (15 cm)	12" Length (30 cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	60" Length (152cm)
	<b>FL6006P</b>	<b>FL6012P</b>	<b>FL6024P</b>	<b>FL6036P</b>	<b>FL6048P</b>	<b>FL6060P</b>

							
Male to Female	6" Length (15 cm)	12" Length (30 cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	48" Length (122cm)	60" Length (152cm)
	<b>FL6106P</b>	<b>FL6112P</b>	<b>FL6124P</b>	<b>FL6136P</b>	<b>FL6148P</b>	<b>FL6148CP*</b>	<b>FL6160P</b>

\*unique coil